



Select “Click here to login” link at bottom. Log on using guest as the user name and region8esc as the password.



To view only, enter recording ID. Choose media type and speed.



Once logged on as a guest, choose “View recording list” link.

Codian IP VCR - Recording list - Microsoft Internet Explorer

Address: http://10.44.36.248/recording_list.html?folder=1

Codian VCR 2220
 host: Region8ESC_VCR login: guest

Home Recordings Help

Home > Top-level folder

Recording list

Top-level folder

Recording ID: <not set>
 Recording console ID: <not set>
 Auto attendant ID: <not set>

Sub folders

| Name |
|-----------------------------------|
| PJC Classes |
| Region 8 Sessions |
| TETN Sessions |
| Texarkana Classes |

No stored recordings

| Name | Status | Numeric ID | H.323 gatekeeper | Length | Recorded at |
|----------------------|--------|------------|------------------|--------|-------------|
| No stored recordings | | | | | |

Page 1 1 2 3 4 5 6 7 8 9 10

Page 1 1 2 3 4 5 6 7 8 9 10

(c) Copyright Codian 2003-2006, [License Information](#)

Choose the sub folder TETN sessions.

Codian IP VCR - Recording list - Microsoft Internet Explorer

Address: http://10.44.36.248/recording_list.html#folder=4

VCR 2220
host: Region8ESC_VCR login: guest

Home Recordings Help

Home > Top-level folder > TETN Sessions

Recording list

Folder "/TETN Sessions"

Recording ID: 8386732
Recording console ID: <not set>
Auto attendant ID: 8386411

Sub folders

Name
<move up>

5 stored recordings

| Name | Status | Numeric ID | H.323 gatekeeper | Length | Recorded at |
|---|--------|------------|------------------|-----------------------|--------------------|
| Watch AI Issues for Dia_26Oct06_14.59 | Idle | | n/a | 3 hours, 12 minutes | Oct 26 2006, 14:59 |
| Watch TETN Building the_02Nov06_18.03 | Idle | | n/a | 3 hours, 0 minutes | Nov 2 2006, 18:03 |
| Watch TETN CPS Update a_16Oct06_18.42 | Idle | 20692 | 20692 registered | 2 hours, 31 minutes | Oct 16 2006, 18:42 |
| Watch TETN Dia de los M_03Nov06_07.48 | Idle | | n/a | 7 hours, 35 minutes | Nov 3 2006, 07:48 |
| Watch TETN E. Coli News_13Oct06_14.53 | Idle | 22394 | 22394 registered | 39 minutes, 0 seconds | Oct 13 2006, 14:53 |

(c) Copyright Codian 2003-2006, [License Information](#)

Click on the session link name you wish to download.

Codian IP VCR - Recording "AI Issues for Dia_26Oct06_14.59" - Microsoft Internet Explorer

Address: http://110.44.36.248/recording.html?name=AI%20Issues%20for%20Dia_26Oct06_14.59&Folder=4

VCR 2220
host: Region8ESC_VCR login: guest

Home > [Top-level folder](#) > [TETN Sessions](#) > Recording "AI Issues for Dia_26Oct06_14.59"

Recording

Recording "AI Issues for Dia_26Oct06_14.59"

Status

Time of recording Oct 26 2006, 14:59
Duration 3 hours, 12 minutes
Recorded media Streaming H.263 video (352 x 288 pixels, 128k / 384k) with G.711mu audio
H.263 video (560 x 320 pixels, 384k) with Polycom(R) Siren14(TM) audio
Complete recording file size 1.14Gb (1138877148 bytes) [download recording](#)
Projected MPEG1 file size 646Mb [download MPEG file](#)
Gatekeeper state n/a
Active streaming playbacks 0 (0)
Active H.323 playbacks 0 (0)
Active downloads 0 (1)
Completed playbacks 0
Completed downloads 2

[Reset maximum values](#)

(c) Copyright Codian 2003-2006, [License Information](#)

Right click on the download MPEG file link and choose save target as. Choose a folder you wish to save the file in. An example would be to create a folder in my documents and name it TETN.

